

# Wide Wheel Roll Stand VGA-60M

## SPECIFICATIONS

### Fixed Height Roll Stand

- The carrier unit has 4 arms and is made of injection aluminum.
- Ø125 plastic wheels that has high carrying capacity, durable are used.
- Two of the plastic wheels are braked and lockable.
- Stand body is made of extruded aluminum profile.
- The stand body has a carrying capacity of 60kg.
- The dimensions of the wire basket are 320x195x125mm.
- The basket has cable hangerson its front and sides.
- The basket has a carrying capacity of 8kg.
- The stand set is made of a lean metal alloy.
- The stand set is fully compatible with MX700, MX750, MX800, MX850.
- The handle is fixed on the part.
- The handle is coated with plastic.
- There is a module apparatus for module.



## BOX CONTENTS

Subpart	Definition	Box No	QTY
PR.AL.M09	Body Profile	1	1
HM.TS.001	Wire Basket	2	1
MG.78	Stand Set	3	1
EN.AL.M02	Carrier Unit	4	1
KZ.M78	Mounting Plate	3	1

DIEMENSIONS

1/1



