

# STAR WHEEL NST CART

## SMS-30Y-NST

### SPECIFICATIONS

Fixed Height Roll Stand


- 60 cm base,
- 75 mm casters, two casters can be braked,
- Built-in counterweight,
- 87 cm post,
- 320 mm x 195 mm x 125 mm utility basket with cord hooks,
- Handle,
- Overall height - floor to mounting plate is : 100 cm,
- Weight: 11 kg,
- Max Carrying Capacity: 45kg.



### BOX CONTENTS

Subpart	Definition	Box No	QTY
PR.AL.008	Body Profile	1	1
HM.TS.001	Wire Basket	2	1
SC.04.NST	Stand	3	1
EN.AL.001	Star Wheel	4	1
EN.PL.013	Handle	2	1

You can scan the QR code to watch the installation video.



DIEMENSIONS

1/1

